Racestyptine

Solution de rétraction gingivale

Composition

Chlorure d'aluminium hexahydraté 25 % m/V, oxyquinol, excipient hydroalcoolique.

Propriétés

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Solution for gingival retraction

Composition

Hexahydrated aluminium chloride 25 % m/V, oxyquinol, hydroalcoholic excipient.

Properties

Aluminium chloride has an astringent effect. It induces a local haemostasis and a contraction of the upper strata of the free gingiva by precipitating tissular and serum proteins.

Indication

- gingival exudation and gingival haemorrhages,
- gingival retraction by in situ impregnation of a cord.

Instructions for use

- Gingival exudation and gingival haemorrhage

Instil one or two drops of solution onto the gingiva using a pipette or a cotton swab. Wait 1 to 2 minutes and rinse carefully the application site.

- Gingival retraction

First insert a gingival retraction cord (non medical) into the sulcus, then impregnate the cord in situ very carefully with a few drops of the solution using a pipette. Leave the cord in place 3 to 5 minutes, then remove it making sure that no fragment remains into the gingival sulcus then wash carefully.

Contra-indications

Possible allergy to aluminium.

Any history of hypersensitivity or photosensitivity to quinoleic ingredients.

Warnings and precautions while using this medicine

- Instill only the number of drops necessary to obtain the required result.
- Limit, as much as possible, the contact time of the solution with the mucosa surface.
- As soon as the local haemostasis has been obtained or the impression has been performed, wash thoroughly the application area so as not to leave any aluminium salt remains potentially irritating.
- Not to be used for over impregnating a medicinal cord.
- Do not ingest.

RACESTYPTINE solution is irritating to eyes and skin.

RACESTYPTINE is a flammable solution.

Presentations

Bottle of 13 ml + one applicator cap.

Bottle of 45 ml.

For professional use only In the practice of dentistry

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